



Professional Daily Care Instructions for Miracle Dressing™

The Miracle Dressing™ System is designed to remain on the wound and is designed to fall off on its own in 3-5 weeks. **Do not remove or change the dressing.**

Dressing should never be removed unless the doctor in charge dictates the removal.

Soak the dressing at least twice a day with a gauze soaked with tap or bottled water.



IMPORTANT: DO NOT REMOVE DRESSING

TO VIEW AN INFORMATIONAL VIDEO REGARDING APPLICATION AND
CLEANING OF SOILED DRESSINGS VISIT OUR WEBSITE:

WWW.MBETHEALTH.COM/APPLICATION-CLEAN

OR SCAN THIS QR CODE:



Post This Card at Nurses Station

For Assistance or to Re-Order:

Extended Kit 6-Pack: MD06-TNO24

Sacral Kit 6-Pack: MD06-TNO04

(800) 431-6238 • 5:30AM - 5PM PST

www.MBETHealth.com